


## Simple Fiber Optic Cable Retractor



## Simple Fiber Optic Cable Retractor



Fiber optic retractors use fiber optic cables to transmit light from a source outside the body to the tip of the retractor. This light illuminates the surgical field, making it ...



Anthony Products stainless steel and silicon molded retractors & instruments are used in surgery for creating and maintaining working space required for the surgical operation in question.



Get the best deals for Fiber optic Retractors at eBay . We have a great online selection at the lowest prices with Fast & Free shipping on many items!




12 Options As low as \$976.89 View Details Fiber Optic Cable for storz-compatible & male acmi \$1,246.39 View Details




The SOAR System is a comprehensive, easy-to-use tubular retraction option that provides simplified, versatile posterior access to the spine. The System offers excellent stability and brilliant illumination ...



This PDF features our “Blue Monarch” 5mm Rigid One-Piece Instrument (Style A) and our 5mm Rotatable One-Piece Instrument (Style B). Both have flush ports for easy cleaning and are available ...



TEBBETS Style Retractors - Serrated tip for use with most fiber-optic cables ... Add to Enquiry SKU: HS-01001 Category: Plastic Surgery



The LiteClip™ surgical retractor connects directly to your fiber optic cable, delivering intense, cool white light right at the point of surgery—no external light source repositioning required.



Find your retractor with fiber optic easily amongst the 9 products from the leading brands (Medicta Instruments, PRECISION SPINE, ...) on MedicalExpo, the medical equipment specialist for your ...



- This malleable fiberoptic light bundle can illuminate the many difficult-to-access surgical sites in the abdomen and pelvis.
- The captive clamp will attach to many standard table mounted and self ...



Our most popular fiber-optic lighted retractors, Tebbetts retractors are designed to enhance visualization during a wide range of procedures. Low-profile oval shaped light bundle provides optimal ...



Device Description: Fiber-optic cables and Fiber-optic Retractors are designed to deliver maximum light when coupled to a medical grade fiber-optic light source. Our fiber-optic cables can be used with ...

## Contact Us

For more information, pricing, or custom network solutions, please contact us:

Website: <https://hashherbcafe.co.za>

Email: [hello@hashherbcafe.co.za](mailto:hello@hashherbcafe.co.za)

Phone: +27 63 814 7295

Address: 15 Galaxy Road, Linbro Business Park, Johannesburg, 2065, South Africa

This document is for informational purposes only. Specifications subject to change without notice.

